

SafeGuard™ Power System

The unmatched technology of the SafeGuard Power System in the LIFEPAK CR® Plus AED was designed to make device maintenance efficient, reliable and affordable.




Cutaway view of the LIFEPAK CR Plus AED

How it Works

Many AEDs are powered by external batteries that lose capacity over time. The SafeGuard Power System offers a dual layer of security inside the *CR Plus* as the internal battery is kept to its optimal power level via the CHARGE-PAK™ battery charger.

Benefits of the SafeGuard Power System

- The *CR Plus* comes with an internal, rechargeable battery that lasts the life of the product.
- The CHARGE-PAK battery charger and electrode replacement cycles are synchronized to expire on the same date, minimizing maintenance management.
- The CHARGE-PAK battery charger and electrode replacement kit are the lowest cost of any on the market.



For further information, please contact Physio-Control at 800.442.1142 (U.S.), 888.879.0977 (Canada) or visit our website at www.physio-control.com.



Physio-Control Headquarters
11811 Willows Road NE
Redmond, WA 98052
www.physio-control.com

Customer Support
P. O. Box 97006
Redmond, WA 98073
Toll Free 800 442 1142
Fax 800 426 8049

Physio-Control Canada
Medtronic of Canada Ltd
6733 Kitimat Road
Mississauga, ON
L5N 1W3
Tel 888 879 0977
Fax 866 430 6115



Physio-Control, Inc., 11811 Willows Road NE, Redmond, WA 98052 USA
